

**CATALOGO DE PRODUCTOS
LABORATORIO
DETECCION DE ENDOTOXINA**





What makes LAL so important?

The LAL test is the most sensitive, accurate and cost effective test on the market today to detect contaminating endotoxins. This test was first licensed by the FDA in the 1970's, and is now the gold standard. It can detect endotoxin in the parts per billion. That's like finding a grain of sand in an Olympic swimming pool. Prior to LAL, rabbits were used to test for endotoxin by injecting the rabbit with samples of the product being manufactured and waiting two or three days to see if the rabbit developed a fever. Hundreds of thousands of rabbits were required to be held and utilized this way. LAL based assays replace this test with one that is more humane, more accurate, cost effective and can give results in a test tube, in about an hour. There are very few people you are likely to meet in your lifetime who have not benefited from a bacterial endotoxins test.

How does the process of bleeding crabs work?

The process is very similar to when people donate blood. The crabs are checked for good health, placed in a very clean laboratory, where we disinfect a portion of the shell, and carefully insert a sterile needle. The crabs have a sinus in the dorsal aspect of their body just under the shell that holds excess blood, we collect from that region. The way the crabs are held, limits the blood that can be harvested to the dorsal sinus, the majority of the blood which is in the gill area is untouched. Studies have shown that the crabs tolerate this process very well and the overwhelming majority survives.

What types of things are tested with the blood?

The FDA has mandated (it is the law) that all injectable or indwelling materials must be tested for endotoxin contamination before being released for sale. This is to protect the public from products that are not sufficiently free of materials that can make a patient ill from exposure to gram-negative cell wall material. If endotoxin enters your blood stream it can make you sick and possibly even kill you. So the test we manufacture is used for medical devices, such as knee replacements, stents, heart valves, intravenous solutions; and drugs and vaccines like childhood immunizations, insulin, flu vaccine and chemotherapy drugs to name a few. Anything injected or implanted into the human body must be free of endotoxin.

Where do the crabs you bleed come from?

Most of the crabs that come to our facility are from Massachusetts waters, Vineyard Sound, Nantucket Sound, and Buzzards Bay. Fishermen catch them a number of different ways but must follow strict regulations on size, number of crabs harvested, and quota management.

DETECCION DE ENDOTOXINA

METODO KINETIC CROMOGENIC



BioTek ELx808 Absorbance Reader

CATALOGO # 25-315S

The ELx808™ Incubating Absorbance Microplate Reader has been validated as part of Lonza's quantitative endotoxin detection system interfacing with Lonza's WinKQCL™ Endotoxin Detection and Analysis Software.



SAMPLE VIAL, PYROGEN-FREE 100/KIT

CATALOGO # 80-507U

10 ML GLASS DEPYROGENATED SAMPLE CONTAINERS WITH SCREW CAPS TO BE USED IN APPLICATIONS REQUIRING THE ENDOTOXIN TEST. 100/BOX



SAMPLE VIAL, PYROGEN-FREE 25/KIT

CATALOGO # 80-507L

10 ML DEPYROGENATED GLASS SAMPLE CONTAINERS WITH SCREW CAPS TO BE USED IN APPLICATIONS REQUIRING THE ENDOTOXIN TEST. 25/BOX.



Endotoxin Challenge Vials 25 vials/kit

CATALOGO # 192568

ENDOTOXIN CHALLENGE VIALSTM, >1,000 EU/VIAL, 25 VIALS

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

METODO KINETIC CROMOGENIC



Annual WinKQCL™ Support for Cat. No. 25-611

WINKQCL Software Support Version 6

CATALOGO #: 25 611SUP

Annual WinKQCL™ Software
Advanced Support for Cat. No.
25-611



Kinetic-QCL 192 Test Kit

CATALOGO # 50-650U

Kinetic Chromogenic LAL Assay
for endotoxin detection for 192
tests, Assay sensitivity: 0.005 -
50 EU/mL



USB Wired 2D Barcode Scanner

COTALOGO # 25-507

Barcode Scanner



PyroGene 192 Test Kit

CATALOGO # 50-658 U

Recombinant Factor C Endpoint
Fluorescent Assay for endotoxin
testing up to 192 tests.


COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

METODO KINETIC CROMOGENIC



LAL Reagent Water- 30ml
CATALOGO #: W50-640



Pyrogen-Free Test Tubes 10x75mm w/Caps
COTALOGO # N201
Pyrogen-free borosilicate glass test tubes with screw-caps. Recommended for gel-clot testing.



Pyrogen-Free Dilution 13x100mm w/oCaps
CATALOGO # N207
Pyrogen-free borosilicate glass dilution tubes to be used in endotoxin testing applications.



PYROGENT-5000 100 Test Kit
CATALOGO # N383
Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 100 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

METODO KINETIC CROMOGENIC



Annual WinKQCL™ Support for Cat. No. 25-611



WINKQCL Software Support Version 6

CATALOGO #: 25 611SUP

Annual WinKQCL™ Software
Advanced Support for Cat. No.
25-611

USB Wired 2D Barcode Scanner

COTALOGO # 25-507

Barcode Scanner



SAMPLE VIAL, PYROGEN-FREE 100/KIT

CATALOGO #: 25 80-507U

10 ML GLASS DEPYROGENATED SAMPLE
CONTAINERS WITH SCREW CAPS TO BE
USED IN APPLICATIONS REQUIRING THE
ENDOTOXIN TEST. 100/BOX



Endotoxin Challenge Vials 25 vials/kit

1ENDOTOXIN CHALLENGE VIALSTM,
>1,000 EU/VIAL, 25 VIALS

CATALOGO #: 192568

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA


LAL ENDOTOXIN DETECTION



Kinetic-QCL Bulk Kit 2400 Test Kit

CATALOGO #: 50-650H

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL



LAL REAGENT GRADE MULTI-WELL PLATES

CATALOGO # 25-340

96-well individually wrapped microplate specific for endotoxin testing.



PYROGENT Plus 200 Tests 0.125 EU/ml Sens (Gel Clot)

CATALOGO #: N294-125

Gel clot LAL assay for 200 tests. Kit contains endotoxin. Sensitivity: 0.125 EU/ml.



Pyrospense Kit Dispensing Reagent

CATALOGO #: N188

Pyrospense™ Reagent, 5 × 5 mL vials

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



Pyrogen-Free Test Tubes 10x75mm w/oCap

CATALOGO #: N205

Pyrogen-free borosilicate glass test tubes without screw-caps. Recommended for gel-clot testing.



LAL Reagent Water - 100ml

COTALOGO # W50-100

LAL Reagent Water, 100 mL



PYROGENT Plus 64 Tests 0.125 EU/ml Sens (Gel Clot)

CATALOGO #: N283-125

Gel clot LAL assay for 64 tests. Kit contains endotoxin. Sensitivity: 0.125 EU/ml.



Endotoxin, E. coli O55:B5

CATALOGO #: 193783

Endotoxin, 5 vials, intended for use in endotoxin removal system challenges and other spiking studies

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA


LAL ENDOTOXIN DETECTION



Tris Buffer, 50mM Solution For Use w/ LAL

CATALOGO #: S50-642

50 mM Tris Buffer, 30 mL



Reagent Reservoirs, 10/pkg

COTALOGO # 190035

LAL Reagent Reservoirs used with multi-channel pipettors ideal for endotoxin testing.
10/pack.



PYROGENT Plus 200 Tests 0.03 EU/ml Sens (Gel Clot)

CATALOGO #: N294-03

Gel clot LAL assay for 200 tests.
Kit contains endotoxin. Sensitivity:
0.03 EU/ml.



Beta-G-Blocker Kit

CATALOGO #: N190

b-G-Blocker™, 5 x 5 mL

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



Lonza PyroWave XM Reader
CATALOGO #: 25-345S

The PyroWave™ XM Fluorescence Reader brings a new generation in fluorescence technology to users of the PyroGene™ Recombinant Factor C Endotoxin Detection Assay. This fluorescence microplate reader is optimized for Lonza's PyroGene™ Assay, bringing new and improved technology to the laboratory.



PYROGENT Plus 200 Tests O.06 EU/ml Sens
COTALOGO # N294-O6

Gel clot LAL assay for 200 tests. Kit contains endotoxin. Sensitivity: 0.06 EU/ml.



Kinetic-QCL Bulk Kit 2040 Test Kit

CATALOGO #: 50-65ONV

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2040 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL



Endotoxin, E.coli O55:B5
CATALOGO #: N186

Control Standard Endotoxin for PYROGENT™ Gel Clot LAL Assays, E coli Strain O55:B5, 5 vials

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



PYROGENT Plus 64 Tests 0.06 EU/ml Sens
CATALOGO #: N283-06

Gel clot LAL assay for 64 tests. Kit contains endotoxin. Sensitivity: 0.06 EU/ml.



PPYROGENT Plus 200 Tests 0.25 EU/ml Sens
COTALOGO # N284-25

Gel clot LAL assay for 200 tests. Kit contains endotoxin. Sensitivity: 0.25 EU/ml.



PYROGENT-5000 200 Test Kit
CATALOGO #: N384

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 200 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml



LAL Reagent Water, 1000ml
CATALOGO #: N186

LAL Reagent Water, 1000 mL

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com


DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION




Pyrogene Bulk Kit 2880 Test Kit
CATALOGO #: 50-658NV

Recombinant Factor C Endpoint
Fluorescent Assay for endotoxin testing up
to 2,880 tests.



MgCl₂, 10mM Solution For Use with LAL
COTALOGO # S50-641

10 mM MgCl₂ Solution, 30 mL



Pyrogen-5000 Bulk Kit 2250 Test Kit
CATALOGO #: N588

Bulk turbidimetric kinetic LAL Assay for
endotoxin testing for 2,250 tests. Kit
contains endotoxin and reconstitution
buffer. Assay sensitivity: 0.01 to 100
EU/ml.



PYROGENT 80 Tests 0.125 EU/ml
Sens

Gel clot LAL assay for 80 tests. Kit does not
contain endotoxin. Sensitivity: 0.125 EU/ml.

CATALOGO #: N186

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA


LAL ENDOTOXIN DETECTION



EPPENDORF BIOPUR TIPS 1000 UL

CATALOGO #: 25-417

Pyrogen-free Eppendorf® Biopur Tips, 50-1000 µL. Ideal for use in the endotoxin test.



Pyrogen-5000 Bulk Kit 4500 Test KitL

COTALOGO # N688

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 4,500 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.



PPyrogen-5000 Bulk Kit 25 Vials/Kit

CATALOGO #: T50-300L

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 1,250 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.



LAL Reagent Water - 500ml Sens

LAL Reagent Water, 500 mL

CATALOGO #: W50-500

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



EPYROGENT 250 Tests 0.06 EU/ml Sens
CATALOGO #: N194-O6
Gel clot LAL assay for 250 tests. Kit does not contain endotoxin. Sensitivity: 0.06 EU/ml.



PPYROGENT Plus 4000 Tests 0.125 EU/ml Sen
COTALOGO # N494-125
Gel clot LAL assay for 4,000 tests. Kit contains endotoxin. Sensitivity: 0.125 EU/ml.



PYROGENT Plus 4000 Tests 0.06 EU/ml Sens
CATALOGO #: N494-O6
Gel clot LAL assay for 4,000 tests. Kit contains endotoxin. Sensitivity: 0.06 EU/ml.



PYROGENT Plus 4000 Tests 0.03 EU/ml Sens
CATALOGO #: N494-O3
Gel clot LAL assay for 4,000 tests. Kit contains endotoxin. Sensitivity: 0.03 EU/ml.

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



PYROGENT 250 Tests 0.25 EU/ml Sens
CATALOGO #: N184-25

Gel clot LAL assay for 250 tests. Kit does not contain endotoxin. Sensitivity: 0.25 EU/ml.



BioTek Absorbance Test Plate

COTALOGO # 7260522

BioTek™ Absorbance Test Plate



Pyrogen-5000 Bulk Kit 100 Vials/Kit
CATALOGO #: T50-300U

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 5,000 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/m



PYROGENT 80 Tests 0.06 EU/ml Sensitivity

CATALOGO #: N183-06

Gel Clot LAL Assay for endotoxin testing with sensitivity of 0.06 EU/ml, 80 tests. Kit does not include endotoxin.

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



Pyrogen-5000 Bulk Kit 100 Vials/Kit
CATALOGO #: T50-600U

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 10,000 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.



Pyrogen-5000 Bulk Kit 25 Vials/Kit
CATALOGO #: T50-600L

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 2,500 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.



PYROGENT 250 Tests 0.03 EU/ml Sens
COTALOGO # N194-03

Gel clot LAL assay for 250 tests. Kit does not contain endotoxin. Sensitivity: 0.03 EU/ml.



Stepped Neutral Density Test Plate

CATALOGO #: 25-342
Stepped Neutral Density Plate

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

LAL ENDOTOXIN DETECTION



EPPENDORF BIOPUR TIPS 300 UL†

CATALOGO #: 25-416

Eppendorf® Biopur Tips, 20-300 µL, 5x96 tips/tray. Ideal for use in the endotoxin test with a multi-channel pipettor.



PYROGENT Plus 4000 Tests O.25 EU/ml Sens
COTALOGO # N288-25

Gel clot LAL assay for 4,000 tests. Kit contains endotoxin. Sensitivity: 0.25 EU/ml.



WinKQCL Reader License Version 6

CATALOGO #: 25-613

WinKQCLTM 6 Software Reader License



WinKQCL Software Package Version 6

CATALOGO #: 25-611

WinKQCLTM 6 Endotoxin Detection and Analysis Software Package is a fully-integrated, 21CFR part 11 compliant solution for quantitative endotoxin detection analysis, data management, reporting, and trending. (disc copy)

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA


LAL ENDOTOXIN DETECTION




Pyrogen-5000 Lysate Vial 10.4 ml
CATALOGO #: T50-600



PYROGENT 250 Tests 0.125 EU/ml Sens
COTALOGO # N194-125



EPPENDORF BIOPUR TIPS 200 UL
CATALOGO #: 25-415



**USP REFERENCE STANDARD
ENDOTOXIN**
CATALOGO #: E700

COMPROMETIDOS CON LA
CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

Endotoxin Detections kits

DETECCION DE ENDOTOXINA

KITS SENTITIVITY 0.0005 - 50 EU/ML



Kinetic-QCL™ Kinetic Chromogenic LAL Assay, 192 Test Kit

Kinetic Chromogenic LAL Assay for endotoxin detection for 192 tests, Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650H



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 600 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 600 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643L



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay; 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2040 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2040 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

KITS SENTITIVITY 0.005 - 50 EU/ML



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay; 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643U



Kinetic-QCL™ Kinetic Chromogenic LAL Assay, 192 Test Kit

Kinetic Chromogenic LAL Assay for endotoxin detection for 192 tests, Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2040 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2040 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650NV



PyroGene™ Recombinant Factor C Endpoint Fluorescent Assay

Recombinant Factor C Endpoint Fluorescent Assay for endotoxin testing up to 192 tests.

Catalog #: 50-658U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 600 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 600 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643L



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL



PyroGene™ Bulk Recombinant Factor C Endpoint Fluorescent Assay

Recombinant Factor C Endpoint Fluorescent Assay for endotoxin testing up to 2,880 tests.

Catalog #: 50-658NV

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA



829-988-8762



mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA


KITS SENTITIVITY 0.01 - 100 EU/ML



Kinetic-QCL™ Kinetic Chromogenic LAL Assay, 192 Test Kit

Kinetic Chromogenic LAL Assay for endotoxin detection for 192 tests, Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay; 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643U



PyroGene™ Recombinant Factor C Endpoint Fluorescent Assay

Recombinant Factor C Endpoint Fluorescent Assay for endotoxin testing up to 192 tests.

Catalog #: 50-658U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2040 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2040 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650NV



PYROGENT™-5000 Kinetic Turbidimetric LAL Assay, 100 Test Kit

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 100 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.


Catalog #: N383



PYROGENT™-5000 Kinetic Turbidimetric LAL Assay


Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 200 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

Catalog #: N384



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 600 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 600 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643L

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA



829-988-8762



mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA


KITS SENTITIVITY 0.01 - 100 EU/ML



PyroGene™ Bulk Recombinant Factor C Endpoint Fluorescent Assay

Recombinant Factor C Endpoint Fluorescent Assay for endotoxin testing up to 2,880 tests.


Catalog #: 50-658NV



PYROGENT™-5000 Bulk Kinetic Turbidimetric LAL Assay

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 5,000 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.


Catalog #: T50-300U



PYROGENT™-5000 Bulk Kinetic Turbidimetric LAL Assay, 2250 Test Kit

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 2,250 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.


Catalog #: N588



PYROGENT™-5000 Bulk Kinetic Turbidimetric LAL Assay, 10000 Test Kit

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 10,000 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

Catalog #: T50-600U



PYROGENT™-5000 Bulk Kinetic Turbidimetric LAL Assay, 4500 Test Kit

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 4,500 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.


Catalog #: N688



PYROGENT™-5000 Bulk Kinetic Turbidimetric LAL Assay

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 2,500 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

Catalog #: T50-600L



PYROGENT™-5000 Bulk Kinetic Turbidimetric LAL Assay

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 1,250 tests. Kit does not contain endotoxin or reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

Catalog #: T50-300L

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com

DETECCION DE ENDOTOXINA

KITS SENTITIVITY 0.02 EU/ML



Kinetic-QCL™ Kinetic Chromogenic LAL Assay, 192 Test Kit

Kinetic Chromogenic LAL Assay for endotoxin detection for 192 tests, Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay; 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643U



PyroGene™ Recombinant Factor C Endpoint Fluorescent Assay

Recombinant Factor C Endpoint Fluorescent Assay for endotoxin testing up to 192 tests.

Catalog #: 50-658U



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2040 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2040 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650NV



PYROGENT™-5000 Kinetic Turbidimetric LAL Assay, 100 Test Kit

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 100 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

Catalog #: N383



PYROGENT™-5000 Kinetic Turbidimetric LAL Assay

Bulk turbidimetric kinetic LAL Assay for endotoxin testing for 200 tests. Kit contains endotoxin and reconstitution buffer. Assay sensitivity: 0.01 to 100 EU/ml.

Catalog #: N384



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 2400 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 2400 tests. Kit contains endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: 50-650H



Kinetic-QCL™ Bulk Kinetic Chromogenic LAL Assay, 600 Test Kit

Bulk Kinetic Chromogenic LAL Assay for endotoxin detection for 600 tests. Kit does not contain endotoxin. Assay sensitivity: 0.005 - 50 EU/mL

Catalog #: K50-643L

COMPROMETIDOS CON LA CALIDAD EN LA ENTREGA
CONTACTANOS Y COTIZA

☎ 829-988-8762

✉ mrodriguez@ppsolutionsrd.com